

# BIOHIT

Innovating for Health

*Special Offer!*  
*eLINE*  
*Starter Kit*  
*up to 40% off*



Autoclavable  
lower parts



## eLINE<sup>®</sup> PCR Starter Kit

eLINE Starter Kit offers an easy start in PCR applications and all general laboratory work. Biohit eLINE electronic pipettes and dispenser speed up your work and guarantee user-independent results and excellent accuracy even in small volumes. Sterile Biohit SafetySpace Filter Tips and Dispenser Tips ensure contamination-free pipetting.

### eLINE<sup>®</sup> - unbeatable accuracy

- Fully electronic operation and tip ejection
- User-independent results
- Twice as fast as a mechanical pipette

The Biohit Cooling Rack keeps your PCR reagents cool during pipetting. The aluminum rack is autoclavable and it holds 12 micro tubes (1,5 and 2 ml).

### Starter Kit Cat. No. 730820 includes:

eLINE 0.2-10 µl single-channel pipette  
eLINE 10-300 µl single-channel pipette  
eLINE Lite electronic positive displacement dispenser

Charging carousel for 4 pipettes

SafetySpace™ 10 µl Filter Tip rack, sterile  
SafetySpace™ 300 µl Filter Tip rack, sterile

Dispenser Tip 0.5 ml (10 pcs), sterile  
Dispenser Tip 1 ml (10 pcs), sterile

Aluminium cooling rack for 12 micro tubes

Cryo pen (black)

Literature

**Special Kit Price \$1799.00**

Valid October 15 - December 31, 2010. Available to any authorized participating Biohit dealer. Cannot be combined with any other special offers. Pricing subject to change without notice.



## SafetySpace™ Filter Tips for safer pipetting and better accuracy

SafetySpace Filter Tips have more space between the sample and the filter than conventional filter tips. With SafetySpace Filter Tips you do not need to worry about the sample touching the filter whatever pipetting technique or type of liquid you're using. The tips are packed in color-coded single tray racks. They are certified DNase, RNase and endotoxin free and are e-beam pre-sterilized.

Vol.	Range	Length	Packaging Type	Qty	Cat. No
0.1 - 10 µl	31.5 mm	Single Tray	Pre-Sterilized	10x96	790011F
0.5 - 20 µl	51 mm	Single Tray	Pre-Sterilized	10x96	790021F
2 - 120 µl	51 mm	Single Tray	Pre-Sterilized	10x96	790101F
5 - 200 µl	54 mm	Single Tray	Pre-Sterilized	10x96	790201F
5 - 300 µl	52.5 mm	Single Tray	Pre-Sterilized	10x96	790301F
50 - 1000 µl	78 mm	Single Tray	Pre-Sterilized	10x96	791001F
50 - 1200 µl	90 mm	Single Tray	Pre-Sterilized	10x96	791211F

They are particularly useful in:

- pipetting foaming liquids like buffers and proteins
- pipetting solvents
- reverse pipetting
- multiple dispensing functions of electronic pipettes



## New Biohit Tip Refill System: ValuePak

The new refill system is designed for today's demanding laboratories and scientists. The improved packaging has been designed to be environmentally-friendly with reusable tip boxes and completely recyclable tip trays. The result is less plastic waste, greater convenience, and more savings. Tip boxes are tongue-and-groove for easy stacking.

- Universal design fits most major brands of pipettes
- High quality design following ISO/CE/IVD
- Fully autoclavable
- All sizes from 0.1 - 1200ul



Cat No.	Volume	Length	Packaging	Type	Qty	List Price	Special Price
790012-OP	0.1 - 10 µl	31.5 mm	Refill Tower	Non Sterile	10x960 Tips	\$432.00	Call
790302-OP	0.5 - 300 µl	51 mm	Refill Tower	Non Sterile	10x960 Tips	\$330.00	Call
790352-OP	5 - 350 µl	54 mm	Refill Tower	Non Sterile	10x960 Tips	\$380.00	Call
791002-OP	10 - 1000 µl	71.5 mm	Refill Pack	Non Sterile	10x960 Tips	\$330.00	Call
791202-OP	50 - 1200 µl	71.5 mm	Refill Pack	Non Sterile	10x960 Tips	\$422.00	Call

